

iPhysio becomes *universal* on its 10th anniversary

LYRA ETK is a French medtech company specializing in dental restoration. It combines traditional implantology and the latest technologies to deliver innovative solutions, working toward biological and digital implantology practice.

iPhysio is a unique solution developed by LYRA ETK as a new standard in modern implantology.

In brief, iPhysio is a one piece component that enables the management of the entire restorative sequence, from implant through to delivery of the final crown, with no removal along the way. It provides a true 3-in-1 solution for healing abutment, scanbody, and temporary crown.

Beyond the basics, iPhysio offers much more in reality:

- Its U-shape base helps prevent gin-

gival crest recession

- Its bioactive antibacterial surface promotes soft tissue healing
- Morphological shapes produce the perfect emergence profile
- As iPhysio celebrates its 10th anniversary, it is becoming universal

iPhysio is now:

- Recognised by all intraoral scanners on the market
- Compatible with major implant brands in the market
- Expanded in range to tackle all clinical situations: 5 shapes, 4 heights, 2 angulations

Reinforcing iPhysio's excellence, it is now endorsed by the Digital Dentistry Society (DDS) as meeting the highest standard of quality and performance. Clinicians can confidently include it in their practice and digital workflow.



FIG. 1 Lyra ETK iPhysio.

From a scientific standpoint, the results of a retrospective clinical study—covering 103 patients over one year—were released in the Journal of Dentistry exactly a year ago: The results were exceptional in terms of bone and soft tissue retention, as well as for biological complications, which featured in less than 1% of cases.

Happy 10th anniversary iPhysio!
Join the new standard of modern implantology.